



TOC TNTplus Vial Test, HR (30-300 mg/L C)

Product #:	TNT811
USD Price:	\$345.00
Available	

Expert High Range TOC test made simple.

TOC TNTplus® chemistry and your Hach spectrophotometer are engineered to simplify water analysis for accurate results, everytime. For determination of Total Organic Carbon (TOC) by the Direct Method.

Easy and safe handling

The innovative vial and reagent delivery makes tests easier to use than powder pillows or liquid reagents and significantly reduces spillage or contamination risk.

No reagent blank necessary

The high quality of TNTplus vials, tight reagent production controls, instrument calibration verification, and high instrument stability all combine to eliminate the need to run reagent blanks—saving you time and money!

Automatic method detection

The bar-coded vials allow for automatic method detection and measurement by the spectrophotometer, considerably increasing the speed of analysis. Errors are reduced as instrument averages 10 readings and rejects outliers from scratched, flawed, or dirty glassware.

Documented shelf life and COA

New 2D barcode details batch number and expiry date of reagents, which are documented along with the measurement result. An automatic warning is issued if expiry date has passed. Lot specific information (COA) is available on RFID tag on the box.

Specifications

Digestion Required:	Yes
EPA compliant:	Yes
Instrument:	DR3900, DR6000, DR1900, DR2800, DR3800, DR5000
Method:	10267
Method Name:	Direct Method
Number of tests:	25
Parameter:	TOC (Total Organic Carbon)
Platform :	TNT plus™
Range:	30 - 300 mg/L C
Shelf Life:	12 months from production date

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)

- TOC X5 Shaker (Item LQV148.99.00002)