MATERIAL SAFETY DATA SHEET



Buffer Solution, pH 9

SECTION 1. Product and Company Idenfication

Product Name and Synonym: Buffer Solution, pH 9

Product Code: 1520

Material Uses:

Manufacturer: Aqua Solutions, Inc

6913 Hwy 225


Deer Park, TX 77536

(281) 479-2569

Entry Date : 12/1/2014
Print Date: 12/1/2014

24 Hour Emergency Assistance : Chemtrec 800-424-9300

Canutec 613-996-6666



SECTION 2 HAZARD IDENTIFICATION

Generally not hazardous in normal handling, however good laboratory practices should always be used. Avoid long term exposure to skin or by inhalation.

SECTION 3 MIXTURE COMPONENTS

SARA 313 C	omponent	CAS Number	Percent Comp.	Dimension	Exposure Limits
☐ Boric Acid		CAS# 10043- 35-3	0.3%	W/V	OSHA PEL 15 mg/mf (total dust)
☐ Potassium (Chloride	CAS# 7447-40- 7	0.4%	W/V	None Established
✓ Sodium Hyd	droxide	CAS# 1310-73- 2	<0.1%	W/V	OSHA PEL 2 mg/mf ACGIH 2mg/mf
☐ Water, Deid Type II	nized ASTM	CAS# 7732-18- 5	Balance	V/V	None Established

SECTION 4 FIRST AID MEASURES

Generally not hazardous in normal handling, however good laboratory practices should always be used. Avoid long term exposure to skin or by inhalation.

FIRST AID: SKIN: Wash exposed area with soap and water. If irritation persists, seek medical attention.

EYES: Wash eyes with plenty of water for at least 15 minutes, lifting lids occasionally. Seek Medical Aid. INHALATION: Remove to fresh air. If not breathing, give artificial respiration. If breathing is difficult, give oxygen

INGESTION: If swallowed, induce vomiting immediately after giving two glasses of water. Never give anything by mouth to an unconscious person.

Buffer Solution, pH 9

SECTION 5 FIRE FIGHTING MEASURES

Any means suitable for extinguishing Fire Extinguisher Type:

None

surrounding fire

Fire / Explosion Hazards:

Wear self-contained breathing apparatus and protective clothing Fire Fighting Procedure:

to prevent contact with skin and clothing.

SECTION 6 ACCIDENTAL RELEASE MEASURES

Absorb spill with inert material. Flush excess down drain.

SECTION 7 HANDLING AND STORAGE

This Material is not considered hazardous. Handle using

safe laboratory practices.

SECTION 8 EXPOSURE CONTROLS/PERSONAL PROTECTION

None required **Respiratory Protection:**

Ventilation

Local Exhaust 🗸

Mechanical

Protective Gloves: Gloves to prevent exposure

Eye Protection: Safety Glasses

w/Side Shields

Use safe laboratory handling Other Protective Equipment:

procedures.

SECTION 9 PHYSICAL AND CHEMICAL PROPERTIES

>95% <5°C Percent Volatile by Volume: **Melting Point:**

>95°C Information not **Evaporation Rate Boiling Point:**

available

not available Vapor Pressure: **Evaporation Standard**

not available Vapor Density: **Auto Ignition Temp** Not applicable Soluble Solubility in Water: Lower Flamm, Limit in Air Not applicable

Clear Appearance /Odors: Upper Flamm. Limit in Air Not applicable

colorless odorless liquid

Flash Point: Not applicable

Specific Gravity:

SECTION 10 STABILITY AND REACTIVITY INFORMATION

Stability: Stable

Conditions to Avoid: Protect from freezing.

Materials to Avoid: Strong Acids, Bases. Water reactive materials

Hazardous Decomposition

Products:

Not known to occur

Hazardous polymerization: Will Not Occur Buffer Solution, pH 9

Conditions to Avoid: None known

SECTION 11 Toxicological Information

Routes of Entry: Absorbed through skin. Eye contact.

Toxicity to Animals:

LD50: [Rat] - Route: oral; Dose: > 90 ml/kg LC50: Not available.

Chronic Effects on Humans: Not available.

Other Toxic Effects on Humans:

Non-corrosive for skin. Non-irritant for skin. Non-sensitizer for skin. Non-permeator by skin.

Non-hazardous in case of

ingestion. Non-hazardous in case of inhalation. Non-irritant for lungs. Non-sensitizer for lungs.

Non-corrosive to the eyes. Noncorrosive

for lungs.

Special Remarks on Toxicity to Animals: Not available. Special Remarks on Chronic Effects

on Humans: Not available.

Special Remarks on other Toxic Effects on Humans: Not available.

SECTION 12 Ecological Information

Routes of Entry: Absorbed through skin. Eye contact.

Toxicity to Animals:

LD50: [Rat] - Route: oral; Dose: > 90 ml/kg LC50: Not available.

Chronic Effects on Humans: Not available.

Other Toxic Effects on Humans:

Non-corrosive for skin. Non-irritant for skin. Non-sensitizer for skin. Non-permeator by skin.

Non-hazardous in case of

ingestion. Non-hazardous in case of inhalation. Non-irritant for lungs. Non-sensitizer for lungs.

Non-corrosive to the eyes. Noncorrosive

for lungs.

Special Remarks on Toxicity to Animals: Not available. Special Remarks on Chronic Effects

on Humans: Not available.

Special Remarks on other Toxic Effects on Humans: Not available.

SECTION 13 Disposal Considerations

Waste must be disposed of in accordance with federal, state and local environmental control

regulations

SECTION 14 Transport Information

DOT Classification: Not Regulated

DOT Regulations may change from time to time. Please consult the most recent D.O.T. regulations.

SECTION 15 Regulatory Information

Federal and State Regulations: TSCA 8(b) inventory: Water

Other Regulations: EINECS: This product is on the European Inventory of Existing

Commercial Chemical Substances.

Other Classifications:


WHMIS (Canada): Not controlled under WHMIS (Canada).

DSCL (EEC):

This product is not classified according to the EU regulations. Not applicable.

SECTION 16 Additional Information

Revisions



3/13/2013	0.1	Revised to 16 sections. LS
3/4/2009	0	Creation date LS

The information herein is believed to be accurate and is offered in good faith for the user's consideration and investigation. No warranty either expressed or implied is made for the completeness or accuracy of the information whether originating from the above mentioned company or not. Users of this material should satisfy themselves by independent investigation of current scientific and medical knowledge that the material can be used safely.