



HQ440d Benchtop Meter Package with ISENA381 Sodium ISE Electrode

Product #: 8507000Call for ship date

Designed for your water applications, the Hach HQ440d laboratory meter is an advanced meter that takes the guesswork out of measurements. HQd^{TM} meters connect with IntelliCALTM probes that automatically recognize the testing parameter, calibration history, and method settings to minimize errors and setup time.

Flexible benchtop meter makes your lab testing of critical water quality parameters more efficient

Two input channels for simultaneous measurement of pH, Conductivity, Dissolved Oxygen, BOD, ORP, Ammonia, Ammonium, Fluoride, Chloride, Sodium, and temperature - any IntelliCAL™ smart probe.

Intuitive user interface for simple operation and accurate results

Guided calibration and check standard routines reduce calibration errors. Stabilization alerts and visual measurement lock ensure that you can trust in the accuracy of the results.

Trust your measurements - IntelliCAL™ smart probes store all calibrations in the probe

Calibration history allows quick and easy change out of probes without re-calibrating. The HQd smart system records serial numbers, current calibration data, user ID, sample ID, time, and date automatically in the data log for complete GLP traceability.

Enhance productivity with an ultra-bright screen and large font size for easy to read results

Internal USB ports simplify data transfer, peripherals, and software updates and connection to PCs printer and keyboards.

Convenient kit includes everything you need to start testing

Meter package includes HQ440d Benchtop Meter, ISENA381 sodium electrode, sodium standard solutions, ionic strength adjustor powder pillows, probe stand, 4 AA batteries, power adapter, USB cable for data transfer, quick-start guide, user manual, and documentation CD

Specifications

Conductivity resolution: $0.01 \mu \text{S/cm}$ with 2 digits

Inputs: 2

IP Rating: IP67

Parameter: pH, mV, ISE, DO, Conductivity, TDS, Salinity, Resistivity, ORP, Temp

What's in the box?

Meter package includes HQ440d Benchtop Meter, ISENA381 sodium electrode, sodium standard solutions, ionic strength adjustor powder pillows, probe stand, 4 AA batteries, power adapter, USB cable for data transfer, quick-start guide, user manual, and documentation CD