



Intellical™ ISECL181 Chloride Ion Selective Electrode (ISE), 1 meter cable

Product #: ISECL18101

IntellicalTM ISECL181 is a digital, combination chloride ion selective electrode (ISE) with a non-refillable, double junction reference and built-in temperature sensor. The Intellical ISECL181 probe is available with a 1 or 3 meter cable and is intended for laboratory use. The ISECL181 is ideal for measuring chloride concentrations in wastewater, drinking water and general water quality applications. The ISECL181's solid-state sensor design eliminates membrane replacement and allows for dry storage of the ISE.

Requires virtually no maintenance

The solid-state sensor design allows for DRY storage of the ISE without a shelf life. Also, the probe does not require the electrolyte to be refilled or membrane replacement.

Fast, stable, and accurate results

Intellical digital probes lock on the result when the measurement is stable, removing the guesswork in having a moving number. Additionally, the probe's large annular reference junction provides maximum surface area between sample and reference electrolyte for fast measurement stability.

Ultimate traceability in measurement history

Stored time and date stamp for each measurement, operator and sample ID, calibration history, parameter, and probe serial number

Can be moved between meters without the need to re-calibrate or re-enter measurement settings

Ideal for a multi-user environment with multiple HQD series Laboratory and Portable meters

Specifications

Accuracy: ± 0.02 mV or 0.05%, whichever is greater

Cable Length: 1 m

Connector: HQD specific

Diameter: 12 mm

Dimensions (D x L): 12 mm x 220 mm

Filling Solution: Non-refillable Dritek gel

ISA Required: 2318069

Junction: Porous PTFE Annular Ring

Length: 220 mm

Material: Sensor Body: Epoxy

Parameter: Chloride
Probe Type: Standard
Reference Type: Ag/AgCl

Resolution: 2, 3, 4 significant digits selectable

Sensor Type: Solid-state crystal membrane

Special Feature: Dry storage & fast response time. No replacement membranes.

What's in the box?

Includes: IntelliCAL ISECL181 combination Chloride ISE with integrated temperature sensor, calibration certificate, and Basic User Manual.

Required Accessories

• Chloride Ionic Strength Adjustor (ISA) Buffer Powder Pillows, pk/100 (Item 2318069)