



# IntelliCAL™ ISEF121 Fluoride Ion Selective Electrode (ISE), 1 meter cable

Product #: ISEF12101

IntelliCAL<sup>TM</sup> ISEF121 is a digital, combination fluoride ion selective electrode (ISE) with a non-refillable, single junction reference and built-in temperature sensor. The IntelliCAL ISEF121 probe is available with a 1 or 3 meter cable and is intended for laboratory use. The ISEF121 is ideal for measuring fluoride concentrations in drinking water, wastewater, and general water quality applications. The ISEF121's solid-state sensor design eliminates membrane replacement and allows for dry storage of the ISE.

## Requires virtually no maintenance

The solid-state sensor design allows for DRY storage of the ISE without a shelf life. Also, the probe does not require the electrolyte to be refilled or membrane replacement.

#### Fast, stable, and accurate results

IntelliCAL digital probes lock on the result when the measurement is stable, removing the guesswork in having a moving number. Additionally, the probe's large annular reference junction provides maximum surface area between sample and reference electrolyte for fast measurement stability.

#### Ultimate traceability in measurement history

Stored time and date stamp for each measurement, operator and sample ID, calibration history, parameter, and probe serial number

#### Can be moved between meters without the need to re-calibrate or re-enter measurement settings

Ideal for a multi-user environment with multiple HQd series Laboratory and Portable meters

### **Specifications**

Accuracy:  $\pm 0.02$ mV or 0.05%, whichever is greater

Cable Length: 1 m

Connector: HQD specific

Diameter: 12 mm

Dimensions (D x L): 12 mm x 220 mm

Filling Solution: Non-refillable driTEK gel

ISA Required: 258999

Junction: Porous PTFE Annular Ring

Length: 220 mm

Material: Sensor Body: Epoxy

Parameter: Fluoride
Probe Type: Standard
Reference Type: Ag/AgCl

Resolution: 2, 3, 4 significant digits selectable
Sensor Type: Solid-state crystal membrane

Special Feature: Dry storage & fast response time. No replacement membranes.

# What's in the box?

Includes: IntelliCAL ISEF121 combination Fluoride ISE with integrated temperature sensor, calibration certificate, and Basic User Manual.