



# IntelliCAL™ Sodium ISE Analysis Package

Product #: ISENA381AP

The IntelliCAL™ Sodium ISE Analysis Package contains the IntelliCAL™ ISENA381 Sodium ISE, 100 pack of sodium ISA powder pillows (4451569), 1 L of 100 mg/L sodium standard solution (2318153), 500 mL of 1000 mg/L sodium standard solution (1474949), 0.02 M ammonium chloride electrode filling solution (2965150), fill solution dispensing kit (2965700).

IntelliCAL™ ISENa381 is a digital, combination sodium ion selective electrode (ISE) with a refillable double junction reference and built-in temperature sensor. The IntelliCAL ISENa381 probe is available with a 1 or 3 meter cable and is intended for laboratory use. The ISENa381 is ideal for measuring sodium concentrations in wastewater, drinking water and general water quality applications.

#### Provides fast, stable, and accurate response

Refillable, easy to maintain design with a unique built-in-temperature sensor. IntelliCAL<sup>TM</sup> digital probes lock on the result when the measurement is stable, removing the guesswork in having a moving number.

#### Ultimate traceability in measurement history

Stored time and date stamp for each measurement, operator and sample ID, calibration history, parameter, and probe serial number.

#### User alerted when re-calibration is needed

No guesswork needed to provide the most reliable and accurate results with each measurement.

#### Can be moved between meters without the need to re-calibrate or re-enter measurement settings

Ideal for a multi-user environment with multiple HQd series Laboratory and Portable meters.

## **Specifications**

Cable Length: 1 m

Diameter: 12 mm

Dimensions (D x L): 12 mm x 220 mm

Filling Solution: 2965126
ISA Required: 4451569
Junction: Ceramic pin

Length: 220 mm

Material: Sensor Body: Zeonor<sup>TM</sup>

Parameter: Sodium
Probe Type: Standard
Sensor Type: Glass

Special Feature: Fast response time

Temperature Range: Continuous use: 0 to 50 °C

## What's in the box?

Kit includes: IntelliCAL ISENa381 combination Sodium ISE (no. ISENA38101), two Hach sodium standard solutions (100 mg/L (1 L, no. 2318153) and 1000 mg/L (500 mL, no. 1474949)), Hach sodium ionic strength adjustor (ISA) Powder Pillows (pack of 100, no. 4451569), 50-mL beaker (no. 108041), and ISE Basic User Manual.

# **Required Accessories**

• Sodium/Potassium Ionic Strength Adjustor Powder Pillows, pk/100 (Item 4451569)