



IntelliCALTM Ammonia ISE Electrode with 1m Cable

Product #: ISENH318101

IntelliCALTM ISENH3181 is a digital, combination gas-sensing ammonia ion selective electrode (ISE) with a replaceable membrane module, refillable outer body, double junction reference and built-in temperature sensor. The IntelliCAL ISENH3181 probe is available with a 1 or 3 meter cable and is intended for laboratory use. The ISENH3181 is ideal for measuring ammonia concentrations in wastewater, drinking water and general water quality applications.

Provides fast, stable, and accurate results

Refillable outer body, an easy to maintain modular membrane design, and a large annular reference junction for maximum surface area between sample and reference electrolyte provide great measurement stability. IntelliCAL digital probes lock on the result when the measurement is stable, removing the guesswork in having a moving number.

Ultimate traceability in measurement history

Stored time and date stamp for each measurement, operator and sample ID, calibration history, parameter, and probe serial number

User alerted when re-calibration is needed

No guesswork needed to provide most reliable and accurate results with each measurement.

Can be moved between meters without the need to re-calibrate or re-enter measurement settings

Ideal for a multi-user environment with multiple HQd series Laboratory and Portable meters

Specifications

Accuracy: ± 0.02 mV or 0.05%, whichever is greater

Cable Length: 1 m

Connector: HQD specific

Diameter: 12 mm

Dimensions (D x L): 12 mm x 220 mm

Filling Solution: 4447226 ISA Required: 4447169

Junction: Porous PTFE Annular Ring

Length: 220 mm

Material: Sensor Body: Epoxy

Parameter: Ammonia
Probe Type: Standard
Reference Type: Ag/AgCl

Resolution: 2, 3, 4 significant digits selectable

Sensor Type: Gas sensing: glass w/ replaceable NH_3 sensitive membrane

Special Feature: Easy-to-replace membrane modules.

Temperature Range: 0 to 50 °C

What's in the box?

Includes: IntelliCAL ISENH3181 combination Ammonia ISE with integrated temperature sensor, 3-pack of replacement IntelliCAL ammonia membrane modules (no. 5812711), electrode fill solution (50 mL, no. 4447226), calibration certificate, and Basic User Manual.

Required Accessories

• IntelliCAL ISENH3181 Ammonia Ion Selective Electrode (ISE) Replacement Membranes, pack of 3 (Item 5812711)